

[Print](#)

## Oral-B SmartSeries PRO 5000

The Oral-B SmartSeries PRO 5000 electric toothbrush with CrossAction toothbrush head is inspired by professional cleaning tools to deliver a superior clean compared to a regular manual toothbrush. Set at a dynamic 16-degree angle for tooth-by-tooth brushing, CrossAction bristles are able to effortlessly attack, lift and sweep away more dental plaque than a regular manual toothbrush.

- Removes up to 100% more plaque than a regular manual toothbrush
- Oral-B 3D movement cleaning system oscillates, rotates and pulsates to break up and remove more plaque than a regular manual toothbrush
- Rechargeable with 5 premium features plus Daily clean: Gum care, Sensitive, Whitening, Deep clean, Wireless SmartGuide and Daily clean
- Senses when you brush too hard with pressure sensor alert
- CrossAction round brush head with perfectly angled bristles for a precise clean
- Bluetooth Connectivity



Smartphone is not included

2 year warranty.

Included in pack: Rechargeable toothbrush handle, Charging station with toothbrush head storage, 3 toothbrush heads: 1 CrossAction, 1 FlossAction, 1 Sensitive, 1 SmartGuide, 1 Travel Case

### Read articles related to this product:

[Power Toothbrush or Manual Toothbrush?](#)

[Using a Rechargeable Electric Toothbrush](#)

[How to Use a Rechargeable Toothbrush. Flossing is not enough. Are you brushing your teeth properly?](#)

### Oral-B SmartSeries PRO 5000 Features:

CROSSACTION® Perfectly angled to get deep in between teeth



#### 5 premium features plus Daily clean:

##### Gum care

Designed to effectively stimulate and massage gums.



##### Sensitive

Operates at a lower speed than Daily Clean to provide a more gentle brushing sensation for sensitive teeth.



##### Whitening

Designed to alternate between speeds for a polishing effect on the tooth surface, which has been demonstrated to be effective at reducing surface teeth stains.



##### Deep Clean

Designed for exceptional mouth cleaning, featuring 45-second quadrant cleaning for total brushing time of 3 minutes, enhancing the results of the standard 30-second quadrant clean of 2 minutes.



##### Wireless SmartGuide

A SmartGuide for personalised Oral Care coaching Displays time, brushing areas and modes.



##### Daily Clean with 3D Movement

Performance mode for daily usage. Oral-B 3D movement cleaning system oscillates, rotates, and pulsates to remove up to 100% more plaque than a manual toothbrush.



##### Brushing Timer

Dentists recommend brushing for two minutes. With the SmartSeries Timer, that goal is more achievable than ever.



##### Brushing Sensor Technology

When you brush too hard, the SmartGuide activates a visual warning, encouraging you to decrease pressure to protect your gums.



##### Battery

A fully charged Oral-B PRO 5000 rechargeable electric toothbrush lasts for up to 10 days with regular use. (when brushing 2 minutes, twice per day).